

# Vitae Mini Lids

## Automatic Capping Machine



### Great Helper

*-Capping diverse sample tubes or centrifuge tubes*



1. **Safe:** can be used in bio-safety hood
2. **Portable:** 24cm×16cm×49cm, 8kg
3. **Convenient:** one hand operation
4. **Efficient:** 4 seconds/sample in average
5. **Intelligent:** automatically trigger lid opening and closing by sensor
6. **Multiple modes:** capping speed and number of rotations can be adjusted as required
7. **Multiple compatibility:** preset 10mL/15mL/50mL and can be adjusted from 0.5mL-50mL



**Raykol Group (XiaMen) Corp., Ltd.**

Add: 5-6F, No.176 Xinfeng Road, HuizhiZone, Torch High-tech Zone, Xiamen, China

Tel:+86-592-5800190

Mail: emma@raykolgroup.com

www.raykolgroup.com